***Nikon***

**Microscope**

**ECLIPSE**

***Ei***

**Instructions**

## Introduction

Thank you for purchasing a Nikon product.

This instruction manual is intended for users of the Nikon Microscope ECLIPSE Ei.

To ensure correct usage, read this manual carefully before operating this product.

- No part of this manual may be reproduced or transmitted in any form without prior written permission from Nikon.
  - The contents of this manual are subject to change without notice.
  - The equipment described in this manual might differ from the actual product in its appearance.
  - Although every effort has been made to ensure the accuracy of this manual, errors or inconsistencies might remain. If you notice any points that are unclear or incorrect, please contact your local Nikon representative.
  - Some of the equipment described in this manual may not be included in the set you have purchased.
  - If you intend to use any other equipment with this product, read the manual for that equipment too.
  - If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment might be impaired.
-

## Contents

Safety Precautions .....	3
Meaning of Symbols Used on the Product.....	3
WARNING and CAUTION Symbols.....	3
⚠ WARNING .....	4
⚠ CAUTION .....	6
Sicherheitsvorkehrungen.....	7
Consignes de sécurité .....	11
Precauciones de seguridad.....	15
Precauzioni per la sicurezza .....	19
Veiligheidsvoorschriften.....	23
Sikkerhedsforanstaltninger .....	27
Precauções de segurança .....	31
Turvatoimet.....	35
Προφυλάξεις Ασφαλείας.....	39
Notes on Handling the Product .....	43
<b>Chapter 1 Names of Parts .....</b>	<b>44</b>
<b>Chapter 2 Observation Steps.....</b>	<b>46</b>
<b>Chapter 3 Detailed Explanation.....</b>	<b>51</b>
3.1 Moving, Installing, and Storing the Microscope.....	51
3.2 Moving the Specimen.....	52
3.3 Focusing.....	53
3.4 Adjustment of Diopter.....	55
3.5 Adjustment of Aperture Diaphragm.....	56
3.6 Oil-Immersion Observation .....	57
3.7 Torque Adjustment of the Coarse Focus Knob .....	58
3.8 Replacing a Specimen Using the Stage Vertical Movement Stopper.....	58
3.9 Optical Characteristics .....	59
3.10 Microscope Terminology .....	59
<b>Chapter 4 Assembly .....</b>	<b>60</b>
<b>Chapter 5 Troubleshooting .....</b>	<b>64</b>
5.1 Optical .....	64
5.2 Operational.....	66
5.3 Electrical.....	67
<b>Chapter 6 Care and Maintenance .....</b>	<b>68</b>
6.1 Cleaning the Lenses .....	68
6.2 Cleaning the Microscope.....	68
6.3 Disinfecting the Microscope .....	68
6.4 Storing the Microscope.....	68
6.5 Periodical Inspections (Charged).....	68

<b>Chapter 7 Specifications .....</b>	<b>69</b>
7.1 Microscopy (Principles) .....	69
7.2 Performance Properties .....	69
7.3 Physical Properties .....	70

## Safety Precautions

Although this product is designed and manufactured to be completely safe during use, incorrect usage or failure to follow the safety instructions provided may cause personal injury or property damage. To ensure correct usage, read this manual carefully before using the product. Do not discard this manual and keep it handy for easy reference.

### Meaning of Symbols Used on the Product

Symbols that appear on the product indicate the need for caution at all times during use. Always refer to the instruction manual and read the relevant instructions before manipulating any part to which the symbol has been affixed.

Symbol	Contents
	<p><b>Biohazard</b></p> <p>This symbol can be found on the microscope body (stage right side), and cautions the following:</p> <ul style="list-style-type: none"> <li>• <b>WARNING:</b> The product may become biohazardous if a sample is spilled onto the product.</li> <li>• To avoid exposure to biohazard, do not touch contaminated parts with your bare hands.</li> <li>• Decontaminate the contaminated parts according to the standard procedures for your facility.</li> </ul>

### WARNING and CAUTION Symbols

Safety instructions in this manual are marked with the following symbols to highlight their importance. For your safety, always follow the instructions marked with these symbols.

Symbol	Contents
 <b>WARNING</b>	Disregarding instructions marked with this symbol may lead to serious injury or death.
 <b>CAUTION</b>	Disregarding instructions marked with this symbol may lead to injury or property damage.

**WARNING****1. Do not disassemble.**

Disassembling this product may result in electric shock or malfunction. Malfunctions and damage due to disassembly will not be warranted. Do not disassemble any part unless instructed to do so in this manual. If you experience problems with the product, contact your nearest Nikon representative.

**2. Read the instructions thoroughly.**

To ensure safety, thoroughly read this manual and the manuals for other equipment to be used with this product. In particular, be sure to follow the warnings and cautions at the beginning of the manuals.

**3. Input voltage, AC adapter, power cord**

The AC adapter for supplying power to this product can be used with 100 to 240 VAC at 50 to 60 Hz. Use only the AC adapter and power cord specified in Chapter 7, "Specifications". Use of another AC adapter or power cord may result in malfunction or fire.

- Avoid using the product in an environment where the supply voltage may fluctuate excessively.
- This product is classified as Class I for electrical shock protection. Be sure to connect it to a protective earth terminal.
- If the AC adapter is covered or items are placed on the AC adapter, heat dissipation may be hindered, causing it to become abnormally hot.

**4. Handling flammable solvents**

The following flammable solvents are used with the product:

- Immersion oil (Nikon Immersion Oil for oil immersion objectives)
- Absolute alcohol (ethyl alcohol or methyl alcohol for cleaning optical parts)
- Petroleum benzine (for wiping off the immersion oil)
- Medical alcohol (for disinfecting the microscope)

Never hold a flame near these solvents. When using a solvent, thoroughly read the instructions provided by the manufacturer, and observe the following precautions.

- Keep solvents away from the product and its surroundings when plugging/unplugging the AC adapter or power cord.
- Be careful not to spill the solvents.

**5. Hazardous samples**

This microscope is mainly for use in microscopic observation of a sample such as cells and tissue affixed to a slide.

When handling a specimen, check to determine whether the specimen is hazardous. Handle hazardous specimens according to the standard procedure for your laboratory. If the specimen is of an infectious nature, wear rubber gloves to avoid infection, and be careful not to touch a specimen. If such a sample comes in contact with the microscope, decontaminate the contaminated portion according to the standard procedure for your laboratory.

## WARNING

### 6. Photobiological safety

This product is designed and manufactured in accordance with the safety standard IEC 62471 “Photobiological safety of lamps and lamp systems”.

Illumination near the stage and light from the camera port of the trinocular eyepiece tube are classified into the following risk group. The distance (hazard distance) from the vicinity of the stage or from the opening of the camera port where the risk group classification is equivalent to the exempt group that does not cause photobiological injury is as follows.

	<b>Risk group classification</b>	<b>Hazard distance</b>
<b>Retinal blue light hazard</b>	Risk group 1	0.5m

Illumination in risk group 1 can be hazardous.

- Pay attention to brightness adjustment by the brightness control knob and avoid looking strong light for a long time.
- Do not look into the field lens. The LED light source is built in under the field lens.
- Do not look into strong light from the camera port from directly above for a long time.
- When you do not attach the camera, attach the supplied cap to the camera port.

**! CAUTION**

**1. Assembly of the microscope**


- Assemble the microscope while the AC adapter is not connected to it.
- Take care to avoid pinching your fingers and hands.

**2. Do not wet the product or allow ingress of foreign matter.**


Do not allow the product to become wet, as doing so may result in malfunction, overheating, or electric shock. If water or other liquids are accidentally spilled onto the product, immediately unplug the power cord from the AC adapter. Then, wipe off the liquid with a dry cloth. Ingress of foreign matter may also result in malfunction. If liquids or foreign matter gets in the product, cease use of the product, and contact your nearest Nikon representative.

**3. Moving the product**

- The microscope has carrying handholds and handles. There are three handholds on the base and two handles on the arm of the microscope. When carrying the microscope, hold the microscope firmly by these handholds and handles. The storage part on the rear of the microscope is not a handle. Do not hold this part as it stores the AC adapter and power cord.
- When moving the product, do not hold by the focus knobs, eyepiece tube, and stage, etc. The parts may become detached and also result in malfunctions.

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Proper way to hold

**4. Do not inadvertently push the microscope on the table.**

If you push this product inadvertently, it might move on the installation surface. We recommend that you attach a wire to the security slot on the back of this product to fix it.

**5. Disposal of the microscope**

To avoid biohazard risks, dispose of the microscope as contaminated equipment, according to the standard procedures for your facility.

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## Notes on Handling the Product

### 1. Handle with care.

The product is a precision optical instrument. Handle the product with care and avoid physical shocks and vibrations.  
In particular, the precision of objectives may be lost by even weak physical shocks.

### 2. Electromagnetic environment

Before using this product, Nikon recommends evaluating the electromagnetic environment of the installation site.

Do not use this product close to strong electromagnetic radiation sources (example: unshielded intentional RF sources.) They may interfere with the proper operation of this product.

This product emits low-level electromagnetic radiation. Do not install this product near precision electronic devices. Otherwise, the performance of such devices might be degraded. If TV or radio reception is affected, move the TV or radio farther away from this product.

### 3. Installation location and storage location

The product is a precision optical instrument. Use or storage of the product under inappropriate conditions may result in malfunctions or loss of precision. Details on the installation and storage environment are given in Chapter 7, "Specifications". Additionally, the following conditions must be considered for the installation location and the storage location.


- Avoid placing the product in a hot or humid location. Use or storage of the product in a hot or humid location may result in molding of or condensation on the lenses, loss of precision, and malfunctions.
- Install the product in a place with little dust and dirt.  
When storing the product, place a cover over the product to protect it from dust.
- Install the product in a place with little vibration.
- Install and store the product on a level and sturdy table or stage that can bear the weight of the product.  
Install the product in a location with minimal exposure to hazards in the event of earthquakes and other potential disasters. If necessary, secure the product to the working desk or other heavy and stable items with a strong wire or other means, so as to prevent it from falling.  
This microscope has a security slot on the back. We recommend that you attach a wire to the security slot to fix the microscope.
- Install the product so that the power cord can be unplugged immediately from the AC adapter in case of an emergency.
- Avoid placing the product in direct sunlight or immediately under room lights.  
The image quality is degraded in a bright environment due to the extraneous light entering the objective. A room light immediately above the microscope may also enter the objective as extraneous light. In this case, it is recommended that you turn off the room light above the microscope.
- Install the product at least 10 cm away from the surrounding walls.
- Do not place items on the product.

### 4. Handling optical parts

Handle optical components such as lenses and filters with care, so as not to damage them. Scratches and dirt (i.e., fingerprints) on the optical components will degrade the microscope image. If components require cleaning, see Chapter 6, "Care and Maintenance".

### 5. Fastening plate


The microscope is held tightly by the red fastening plate during shipment. Be sure to remove the plate before use. The microscope cannot be used without removing the fastening plate. For details, see Section 4.1, "Removing the Fastening Plate".


# 1

## Names of Parts


### Nomenclature of Each Part


### (Fixing part)

### (Holding part, storage part)

**Nomenclature of Control Section**

(1) Diopter ring: Compensates for the difference between the right and left eyesight of the observer.	(12) Status display: Displays the connection state of the power source and the ON/OFF state of the light. Off: The AC adapter (power source) is not connected. Orange: The AC adapter (power source) is connected, and the light ON/OFF switch is off. Green: The AC adapter (power source) is connected, and the light ON/OFF switch is on.
(2) Revolving nosepiece: Switches the objectives.	
(3) Specimen holder: Holds a specimen.	
(4) Stage back & forth movement knob: Moves the specimen holder back and forth.	
(5) Stage right & left movement knob: Moves the specimen holder right and left.	
(6) Stage vertical movement stopper: Sets the upper limit of the stage movement. Used for the fail-safe device for specimen damage protection.	
(7) Adjustment position mark of aperture diaphragm: Display corresponding to the magnification of the objective	(13) Light ON/OFF switch: Turns the light on or off.
(8) Aperture diaphragm lever: Adjusts the aperture diaphragm from $\phi 2$ mm to $\phi 28$ mm.	(14) Brightness control knob: Adjusts illumination. Turned clockwise: Increases brightness. Turned counterclockwise: Decreases brightness.
(9) Fine focus knob: Used for fine focusing.	
(10) Coarse focus knob: Used for coarse focusing.	(15) DC input connector: Connected with the AC adapter.
(11) Coarse focus knob torque adjustment ring: Adjusts the tension of the coarse focus knob.	(16) Security slot: Connected with the security wire.

## 2

## Observation Steps


## 1 Install the microscope.

Carry the microscope from the storage location and place it on a table.


Fix the direction of the eyepiece tube to the observation point (the binoculars are toward you).

- ❗ **Loosen the fixing screw for the eyepiece tube a couple of turns, and then turn the eyepiece tube. → Chapter 4, “Assembly” - Step 2, “Change the direction of the eyepiece tube, remove and attach the eyepiece tube”**

Pull out the AC adapter and power cord from the storage part on the rear of the product.

Change the direction of the eyepiece tube.


Remove the AC adapter.

## 2 Turn on the microscope.

Confirm that the DC plug of the AC adapter is connected to the DC input connector on the rear of the product, and plug the power cord of the AC adapter into the wall outlet.

- ❗ **We recommend that you turn the brightness control knob counterclockwise until it stops before turning on the power. You might feel that the illumination is too bright depending on the position of the aperture diaphragm lever.**


Turn on the microscope.

## 3 Turn on the light.

Check the status.

- Green LED: The light is on.
- Orange LED: The light is off.

If the LED is orange, press the light ON/OFF switch. The LED turns green and the light is turned on.


Turn on the light.

**4 Fully open the aperture diaphragm.**


Set the aperture diaphragm lever to the position of 100x of the adjustment position mark of aperture diaphragm.

**Fully open the aperture diaphragm.**

**5 Adjust the brightness.**

Adjust the brightness by rotating the brightness control knob.



**Adjust the brightness.**

**6 Adjust the interpupillary distance.**

Adjust the distance between the binoculars. (Set the distance so that the right and left viewfields are merged into one when looking through the eyepieces.)

When the binocular part is turned upward, the eye level is raised, and when it is turned downward, the eye level is lowered.


Viewfield

**Adjust the interpupillary distance.**

**7 Align the diopter ring with the standard position.**

Turn the diopter ring on the eyepieces to align its bottom edge with the engraved base line.

Perform this procedure to the right and left eyepieces.


**Set the diopter ring to the standard position.**


## 8 Set the specimen.

Open the claw of the specimen holder with your finger at the root projection. Place the specimen slide on the stage with the coverglass facing upward, and put the claw back to fix it.

Rotate the stage knobs to bring the observation region of the specimen into the optical path. (Let illumination light go through the specimen.) → Chapter 3, “Detailed Explanation” - Section 3.2, “Moving the Specimen”

Set the specimen.


Bring the observation region into the optical path.

## 9 Focus on the specimen with the 10x objective.


Rotate the revolving nosepiece to swing the 10x objective into the optical path. The objective will click into place when rotated into position.

Rotate the coarse and then fine focus knobs to focus on the specimen. → Section 3.3, “Focusing”, and Section 3.7, “Torque Adjustment of the Coarse Focus Knob”


Set the magnification of the objective to 10x.

- ❗ **Do not turn the right and left focus knobs simultaneously in opposite directions. Additionally, do not turn the coarse focus knob further after the stage has reached its upper or lower limit. This operation will damage the focusing mechanism.**

- ✓ **When the focal point cannot be found**  
To find the focal point, rotate the stage vertical movement stopper counterclockwise to unlock the limit.


Focus on the specimen.


### 10 Adjust the diopter.

Swing the 40x objective in the optical path and rotate the focus knobs to focus on the specimen.

Then, change the objective to the 10x (or 4x) objective and focus on the specimen by rotating the diopter ring without using the focus knob. Perform this procedure while looking into the right eyepiece with your right eye and the left eyepiece with your left eye. → Section 3.4, “Adjustment of Diopter”

(1) Focusing with the 40x objective by rotating the focus knobs


(2) Focusing with the 10x objective by rotating the diopter rings

**Adjust the diopter.**

### 11 Change the objective to the magnification used for observation.

Rotate the revolving nosepiece to bring the objective of the desired magnification into the optical path.


**!** Do not directly switch from the 4x objective to the 60x or 100x objective. Be sure to focus on the specimen with the 10x objective before using the objective of a higher magnification.

**Switch to the objective to be used for observation.**

### 12 Adjust the aperture diaphragm.

Set the aperture diaphragm lever to the adjustment position mark of aperture diaphragm (the same figure as the magnification of the objective). → Section 3.5, “Adjustment of Aperture Diaphragm”


**Adjust the aperture diaphragm.**

### 13 Observe the specimen.

When you need to observe the specimen with a different magnification, repeat steps 11 and 12.

### 14 Turn off the microscope.

Set the light ON/OFF switch to OFF (the status display LED turns orange), and unplug the power cord.

**Turn off the microscope.**


## 15 Store the microscope.

After observation, turn the eyepiece tube backward and fix it.

- ❗ **Loosen the fixing screw for the eyepiece tube a couple of turns, and then turn the eyepiece tube. → Chapter 4, “Assembly” - Step 2, “Change the direction of the eyepiece tube, remove and attach the eyepiece tube”**

Store the AC adapter in the storage part on the rear of the microscope, wind the cords around, and then put it back in place.

- ✓ **AC adapter**  
The AC adapter can be stored with it connected to the microscope.

**Fix the eyepiece tube backward and store the AC adapter.**

## 3

## Detailed Explanation


## 3.1 Moving, Installing, and Storing the Microscope

## ■ Moving the microscope

When changing the installation location, or taking the microscope out of the storage location to install it, carry and install it according to the following appropriate method.

The microscope has carrying handholds and handles. There are three handholds on the base and two handles on the arm of the microscope. When carrying the microscope, hold the microscope firmly by these handholds and handles.

- ❗ • The storage part on the rear of the microscope is not a handle. Do not hold this part as it stores the AC adapter and power cord.
- When moving the product, do not hold by the focus knobs, eyepiece tube, and stage, etc. The parts may become detached and also result in malfunctions.


Where you can hold: O  
Where you must not hold: X

## ■ Installing the microscope


Install the microscope in a place appropriate for observation. → “Notes on Handling the Product” - 3, “Installation location and storage location”

When the microscope is carried from the storage location, the binoculars of the eyepiece tube are usually fixed backward. In this case, fix the direction of the eyepiece tube to the observation point (the binoculars are toward you). → Chapter 4, “Assembly” - Step 2, “Change the direction of the eyepiece tube, remove and attach the eyepiece tube”

- ❗ • Loosen the fixing screw for the eyepiece tube a couple of turns, and then turn the eyepiece tube. If the screw head protrudes, turning the eyepiece tube might damage the circular dovetail mount of it by the screw head.

Loosen the screw  
a couple of turns

Pull out the AC adapter and power cord from the storage part on the rear of the microscope.


### ■ Storing the microscope

After using the microscope, arrange it according to the following storage conditions, and then store it in the designated storage location. → “Notes on Handling the Product” - 3, “Installation location and storage location”

Store the AC adapter in the storage part on the rear of the microscope, and wind the cords around. The AC adapter can be stored with it connected to the microscope.

Rotate the eyepiece tube 180 degrees from the observation point and fix it backward.

Then, the upper part of the microscope is made more compact. The binoculars do not hit somewhere, resulting in less storage space. → Chapter 4, “Assembly” - Step 2, “Change the direction of the eyepiece tube, remove and attach the eyepiece tube”

- ❗ **Loosen the fixing screw for the eyepiece tube a couple of turns, and then turn the eyepiece tube. If the screw head protrudes, turning the eyepiece tube might damage the circular dovetail mount of it by the screw head.**


## 3.2 Moving the Specimen

You can set the specimen fixed to the specimen holder to the desired position by operating the stage motion knobs to move the specimen holder back and forth, and right and left. Set the specimen so that the region you want to observe is in the optical path.

The end of the travel scale shows the icon of the stage motion knob to be operated.

### ✓ Travel scale

The stage has the scales that indicate the travel of the specimen in the horizontal and vertical directions. You can use the main scale and vernier scale to read the travel in units of 0.1 mm.



### 3.3 Focusing

This procedure is to operate the fine and coarse focus knobs to raise and lower the stage so that the specimen to be observed can be seen most clearly.

The UP and DOWN indications on the microscope base indicate the direction in which the stage moves by rotation of the focus knobs. For both the fine and coarse focus knobs, turning in the UP direction raises the stage, and turning in the DOWN direction lowers the stage.

#### ✓ Fine scale

One rotation of the fine focus knob moves the stage vertically by 0.2 mm. The fine scale lines are provided at intervals of 2  $\mu\text{m}$  of stage travel.

#### ■ Standard procedure


- (1) Swing the 10x (or 4x) objective into the optical path.
- (2) Rotate the coarse focus knob to raise the stage to its upper limit.  
Since the working distances of 4x and 10x objectives are long, these objectives do not touch the specimen even when the stage is raised to its upper limit provided that the slide and coverglasses of a standard thickness are used. (The standard thickness for slide is 1.2 mm and that for coverglass is 0.17 mm.)
- (3) Looking into the eyepiece, slowly rotate the coarse focus knob to lower the stage. When the specimen image appears, take your hands off the coarse focus knob.

❗ **When rotating the coarse focus knob while looking into the eyepiece, be sure to rotate it in the direction that lowers the stage.**

#### ✓ When the focal point cannot be found

To find the focal point, rotate the movement stopper counterclockwise to unlock the limit.

- (4) Rotate the fine focus knob and precisely focus on the image.

### ■ Using the objective of a higher magnification for observation


First, use the 10x objective to focus on the specimen, and then change it to the desired objective and rotate the fine focus knob to precisely focus on the specimen.

**!** When operating the focus knob while using the objective of a higher magnification, care must be taken to prevent the slide glass from touching the objective. When you raise the stage using the coarse focus knob, look at the gap between the upper surface of the specimen and the front of the objective from the side.

### ■ Using working distances

The stage can be raised by using the working distance (→ See Section 3.10, "Microscope Terminology".) of the objective.

While looking at the microscope from the side, rotate the coarse focus knob with care to bring the specimen close to the objective. Then, when the distance between the specimen and the front of the objective becomes slightly smaller than the working distance, take your hands off the coarse focus knob. The specimen is now almost in focus. Look into the eyepieces and rotate the fine focus knob in the direction that lowers the stage to precisely focus on the specimen.

### ■ Focusing by switching the objective from 40x to 10x

Swing the 40x objective into the optical path. While looking at the microscope from the side, rotate the coarse focus knob until the specimen almost touches the objective (about 0.5 mm apart from the front of the objective). Switch to the 10x objective, look into the eyepieces, and rotate the fine focus knob slightly to find the focal point. Be careful not to hit the objective with the specimen.


### 3.4 Adjustment of Diopter

This procedure is to adjust the diopter ring on the eyepieces according to the difference between your right and left eyesight.


✓ **Effect of diopter correction**

This adjustment enables the user to take full advantage of high performance of objectives, and reduce focus deviation due to switching of the magnification.


- (1) Turn the diopter ring on the right and left eyepieces to align its bottom edge with the engraved base line. This is the standard position for the diopter adjustment.


- (2) Swing the 40x objective in the optical path. Rotate the coarse and then fine focus knobs to bring the specimen in focus. See the description in "■ Focusing by switching the objective from 40x to 10x" in Section 3.3, "Focusing".

- (3) Swing the 10x (or 4x) objective into the optical path.  
 (4) While looking into the right eyepiece with your right eye, focus on the specimen by rotating the right diopter ring without using the focus knob.

- (5) While looking into the left eyepiece with your left eye, focus on the specimen by rotating the left diopter ring without using the focus knob.  
 (6) Repeat steps (2) to (5) to check the specimen is in focus.


### 3.5 Adjustment of Aperture Diaphragm

The size of the condenser aperture diaphragm should be adjusted according to the objective in the optical path.


A guideline for the adjustment position is shown by the “adjustment position mark of aperture diaphragm” on the condenser. Turn the aperture diaphragm lever to “the same figure as the magnification of the objective” in the optical path.

Each time you change the objective magnifications, align the aperture diaphragm lever to the same figure as the magnifying power of the objective.

#### Adjustment position mark of aperture diaphragm

Objective	Color of the mark
4x	Red
10x	Yellow
20x	Light green
40x	Light blue
60x	Dark blue
100x	White/Black

Adjustment position mark of aperture diaphragm

#### ✓ Aperture diaphragm

- The aperture size is increased or decreased by turning the condenser aperture diaphragm lever. If the aperture diaphragm is closed, the brightness and resolution are decreased but the contrast and range of focus are increased. If the aperture diaphragm is opened, the brightness and resolution are increased but the contrast and range of focus are decreased. Since the image resolution will be degraded when the aperture diaphragm is closed too much, do not close the aperture diaphragm too much except when observing a specimen with extremely low contrast, such as a near-transparent specimen.
- The aperture diaphragm controls the numerical aperture of the illumination. Do not use it to control brightness. Use the brightness control knob to control the brightness.

- BE2 PLAN → Objective type
- 40x/0.65 → Magnification 40 times/Numerical aperture 0.65
- ∞/0.17 → Mechanical tube length ∞/Coverglass thickness 0.17 mm
- WD 0.6 → Working distance (Distance between the front of the objective and the specimen) 0.6 mm
- OFN20 → Objective field number 20

#### Specifications of objective

## 3.6 Oil-Immersion Observation

The “Oil” mark on the side of an objective indicates that it is an oil-immersion type objective. (The oil-immersion objective also has a black band around the barrel end.) An oil-immersion objective is used with the immersion oil applied between the front of the objective and the coverglass. For an oil-immersion objective with a numerical aperture of 1.0 or more, also apply immersion oil to the condenser to take full advantage of its performance. (An oil-immersion type condenser needs immersion oil to be applied between the front of the condenser and the specimen.)

The condenser has an oil receptacle around the front of its lens.

### ■ Example of oil-immersion

Apply immersion oil to the condenser and then objective.

#### **Condenser:**

Add a drop of oil on the condenser lens, and then place the specimen on the stage.

#### **Objective:**

Rotate the revolving nosepiece to move the objective out of position. Add a drop of oil to the specimen. Slowly rotate the revolving nosepiece to bring the objective back into position.

### ■ Eliminating air bubbles

Make sure that air bubbles are not trapped during oil application. Air bubbles degrade the image.

Do any of the following to eliminate air bubbles:

- Rotate the revolving nosepiece to move the objective back and forth once or twice.
- Add another drop of oil.
- Wipe off the oil and apply again.

### ■ Handling of the immersion oil

- Use a minimum quantity of oil. If too much oil is applied, surplus oil may flow out onto the stage and the condenser and degrade performance.
- Any oil residue left on the lenses of oil-immersion type objectives or adhesion of oils on the front lens of dry type objectives will degrade image quality. After completing oil-immersion observation, be sure to clean the objective, condenser, and any other parts that may be stained by oil.
- To wipe surplus oil on the front of the condenser lens, loosen the fixing screws for the condenser to remove the condenser.
- Use petroleum benzine to wipe off oil and finish with absolute alcohol (ethyl or methyl alcohol). If petroleum benzine is not available, use methyl alcohol instead. In that case, wipe off the oil several times as the detergency of methyl alcohol is weaker than petroleum benzine. (You can wipe off the oil generally 3 to 4 times.)
- When handling petroleum benzine and absolute alcohol, be sure to follow the instructions described in “4. Handling flammable solvents” in “WARNING” at the beginning of this manual.
- Close the oil container cap tightly after use. Make sure that the cap is closed tight after refilling the container. Check the cap periodically to make sure it has not come loose, allowing oil to leak out.
- Do not press the container hard. Oil may splash out.
- If you find an oil drips around the container, wipe them off.
- Avoid contact of immersion oil with eyes or skin. In the event of contact with eyes or skin, take one of the following measures although Nikon immersion oil does not contain any toxic ingredients.

**Contact with skin: Rinse your skin thoroughly with soap and water.**

**Contact with eye: Rinse your eye thoroughly with water (more than 15 minutes) and see a doctor.**


- Do not leave immersion oil in the sun (ultraviolet rays may damage it).

### 3.7 Torque Adjustment of the Coarse Focus Knob

The tension (torque) of the coarse focus knob rotation can be adjusted.

To increase the tension, turn the coarse focus knob torque adjustment ring counterclockwise.

To decrease the tension, turn the ring clockwise. Do not decrease the tension too much. If it is too loose, the stage will fall under its own weight.

### 3.8 Replacing a Specimen Using the Stage Vertical Movement Stopper


When focusing on the specimen with 40x or higher magnification objective, you will find the specimen is brought very near to the objective<sup>(\*)</sup>.

In such a state, it will be very difficult to replace the specimen while taking care not to break it. In a case like this, use the stage vertical movement stopper for easy specimen replacement.

The stage vertical movement stopper has the indication ("FOCUS LIMIT" "HIGH/LOW") that shows the direction of the rotation and vertical movement. Operate the stage vertical movement stopper according to the following steps.


- (1) Turn the stage vertical movement stopper counterclockwise (in the HIGH direction) until it stops.
- (2) Focus on the specimen, and then turn the stage vertical movement stopper clockwise (in the LOW direction) until it stops. Rotating the stage vertical movement stopper too much causes the set focal point to shift.
- (3) Performing the above steps sets the stage so that it does not move above that position.
- (4) Lower the stage by the focus knob, and replace the specimen. After that, when turning the focus knob, you can find that the stage stops near the position set in (2).
- (5) If the specimen is out of focus, turn the stage vertical movement stopper counterclockwise (in the HIGH direction) to release the set position, and then slightly operate the fine focus knob to adjust the focus.

\*1: The distance between the front of the objective and the specimen when the specimen is in focus is called the "working distance" of the objective. See Section 3.10, "Microscope Terminology".


Turn the stopper until it stops  
(in the HIGH direction)

Turn the stopper until it stops  
(in the HIGH direction)

Focusing

Turn the stopper until it stops  
(in the LOW direction)

Turn the stopper until it stops  
(in the LOW direction)

When raising the stage by the focus knob, the stage stops near the focused position. The stage does not move above that position.

## 3.9 Optical Characteristics

### ■ Combinations of 10x (Field No. 20) eyepiece with BE2-Plan objectives

Objective Magnification [x]	Total Magnification [x]	Numerical Aperture	Real Viewfield [mm]	Depth of Focus [μm]	Resolving Power [μm]	Working Distance [mm]
4	40	0.1	5.0	115.0	2.8	25
10	100	0.25	2.0	18.4	1.1	6.7
20	200	0.4	1.0	6.1	0.7	3.7
40	400	0.65	0.5	2.0	0.4	0.6
60	600	0.8	0.3	1.2	0.3	0.25
100	1000	1.25	0.2	0.7	0.2	0.14

## 3.10 Microscope Terminology

### (1) Total Magnification

The total magnification of a microscope is the individual magnifying power of the objective multiplied by that of the eyepiece.

### (2) Numerical Aperture (N.A.)

The numerical aperture is an important factor in determining the efficiency of the condenser and objective. It is represented by the formula:

$$\text{N.A.} = n \sin \alpha$$

where  $n$  is the refractive index of the medium (air, immersion oil, etc.) between the objective and the specimen or condenser, and  $\alpha$  is the maximum angle at which light enters or leaves the lens from or to a focused object point on the optical axis.

The larger the numerical aperture the brighter the image and the higher the resolution.

### (3) Resolving Power

The ability of an optical system to discriminate between two discrete objects separated by a minute distance. The more minute the distance, the higher the resolving power of the optical system. In relation to the numerical aperture, the resolving power is represented by the following formula:

$$\text{Resolving power} = \frac{\lambda}{2 \times \text{N.A.}} \quad \text{where } \lambda \text{ is the used wavelength of light.}$$

(The resolving power in the above table is indicated for  $\lambda = 0.55 \mu\text{m}$ .)

### (4) Working Distance (W.D.)

The clearance between the front of the objective and the upper surface of the coverglass, when a specimen image is sharply focused. Generally, the higher the magnifying power of the objective, the shorter the working distance.

### (5) Field Number of the Eyepiece

The diameter in mm of the viewfield observable through the eyepiece. When an eyepiece has an indication of "10x / 20", it means that the magnification is 10x and the field number is 20 for that eyepiece.

### (6) Real Viewfield

The diameter of the region of the specimen that is actually observed with the microscope.  
Real viewfield = field number of eyepiece / magnification of objective

### (7) Depth of Focus

The depth (thickness) of the specimen image in focus, extending above and below the focused image plane. The larger the N.A. of the objective, the shallower the depth of focus.

$$\text{DOF} = n \times \left( \frac{\lambda}{2 \times \text{N.A.}^2} + \frac{250000 \times \omega}{M \times \text{N.A.}} \right)$$

DOF: Depth of focus (object side)

$n$ : Refractive index (Varies depending on the medium between the specimen and the objective lens)

For air: 1

For oil immersion (Nikon Immersion Oil): 1.518

$\omega$ : Eye resolution ( $5' = 0.0014$ )

$M$ : Total magnification

$\text{N.A.}$ : Numerical aperture of the objectives

$\lambda$ : Wavelength (observation magnification:  $0.55 \mu\text{m}$ .)

## 4

## Assembly

Before assembling or connecting devices, thoroughly read the "Safety Precautions" and "Notes on Handling the Product" at the beginning of this manual, and be sure to follow all the instructions written therein.

**Tools Required for Assembly:** Hexagonal wrench (the distance between opposite sides is 3 mm)  
(one hexagonal wrench is provided with the microscope)


### 1 Remove the fastening plate

Remove the red fastening plate attached to the top of the stage.

The fastening plate retains movement in the front and back directions of the stage and is attached with two bolts. Remove this fastening plate with the hexagonal wrench provided.

✓ **Keeping the fastening plate**

Be sure to keep the fastening plate with two bolts just in case re-fastening is needed.

### 2 Change the direction of the eyepiece tube, remove and attach the eyepiece tube

When you have purchased the microscope as a set, the eyepiece tube is already installed to the main body.

At the time of observation, turn the binoculars toward you and fix the eyepiece tube. At the time of storage, turn the binoculars backward and fix the eyepiece tube.

When changing the direction of the eyepiece tube, or replacing the eyepiece tube to another one, follow the following procedure.

When installing the eyepiece tube to the main body for the first time, also follow the procedure described in "■ Installing the eyepiece tube".



#### ■ Changing the direction of the eyepiece tube

Manually rotate the fixing screw for the eyepiece tube a couple of turns to loosen it while supporting the eyepiece tube by hand.

Turn the eyepiece tube in the desired direction and securely tighten the fixing screw for the eyepiece tube while the eyepiece tube fits into the circular dovetail mount properly.

✓ **Mark of the direction in which to fix the eyepiece tube for storage**

The microscope has an icon on the arm that indicates the direction in which to fix the eyepiece tube when storing the microscope. This mark indicates that the binoculars of the eyepiece tube is turned backward and fixed.

#### ■ Removing the eyepiece tube

Manually loosen the fixing screw for the eyepiece tube while supporting the eyepiece tube by hand.

Loosen the auxiliary fixing screw for the eyepiece tube by using the supplied hexagonal wrench and remove it.


Remove the eyepiece tube with care not to drop it.

### ■ Installing the eyepiece tube

Insert the eyepiece tube into the circular dovetail mount while tilting it so that the eyepiece tube fits into it properly.

Securely tighten the fixing screw for the eyepiece tube while supporting the eyepiece tube by hand.

Fit the auxiliary fixing screw for the eyepiece tube and tighten it by using the supplied hexagonal wrench.

#### ✓ Function of the auxiliary fixing screw for the eyepiece tube

The auxiliary fixing screw for the eyepiece tube is for preventing the eyepiece tube from falling when the fixing screw for the eyepiece tube is loosened. The eyepiece tube cannot be fastened only with the auxiliary fixing screw for the eyepiece tube.

#### ✓ When screws are attached

If you purchase the microscope without the eyepiece tube, the fixing screw and auxiliary fixing screw for the eyepiece tube are fitted into the screw holes and taped to prevent them from falling. When installing the eyepiece tube, peel off the tape and remove the two attached screws, and then follow the procedure described in “■Installing the eyepiece tube”.

## 3 Attach the camera mount


The camera mount (C mount) is attached to the camera port of the trinocular eyepiece tube.

Loosen the knurled screw on the camera port and remove the cap and C mount. Photomicrography becomes available after attaching the C mount to the camera, and then attaching it to the camera port. The optical path of the trinocular eyepiece tube is fixed at the ratio: binoculars : camera port = 50 : 50.

For details about photomicrography, refer to the instruction manual of the camera.


Do not lose the cap. When you do not attach the camera, make sure that the cap is placed on the camera port.

The cap is attached.

Camera mount (C mount)

Knurled screw


## 4 Install the eyepieces

The 10x (Field No. 20) eyepieces are attached to the eyepiece tube before shipment.

When changing to a 15x (Field No. 18) eyepiece (optional), be sure to change both the right and left eyepieces together. The right and left eyepieces should be of the same magnification.

### ✓ Replacement of eyepieces

To remove the eyepieces, loosen the fixing screws for the eyepieces. The microscope does not come with a hexagonal wrench (the distance between opposite sides is 2.5 mm) to loosen the fixing screws. Contact your nearest Nikon representative if you have any questions.


## 5 Install the objectives

Screw an objective into the revolving nosepiece. Screw the objective firmly to the end.

First, remove the specimen from the stage and lower the stage. Remove the objective while holding it by both hands so that it does not fall, and then install a new objective.

Set the objectives in such an order that the objective magnification increases as the revolving nosepiece is turned clockwise (as viewed from the top of the microscope).

When you have purchased the microscope as a set, three objectives are already installed to the revolving nosepiece. When adding an objective or replacing the objective with another one, also install it in the magnification order as described above.


Example of order of installing the objectives

## 6 Install the condenser

The condenser is attached to the microscope before shipment.


If you remove the condenser and then re-install it, position the condenser with its adjustment position mark of aperture diaphragm facing front. Tighten the two condenser fixing screws.

## 7 Connect the AC adapter and power cord

Plug the power connector of the specified power cord into the inlet for equipment of the AC adapter.

Specified power cord → Chapter 7, "Specifications" - Section 7.3, "Physical Properties"


## 8 Others

### 1) Eyepiece Pointer

The eyepiece pointer (optional) serves as a reference for pointing out the specimen to be observed. Attach the pointer to one of the eyepieces. Rotate and remove the viewfield ring from the eyepiece end. Attach the eyepiece pointer to the eyepiece and then put back the viewfield ring.

### 2) Filter Receptacle and Filters

Two  $\phi 45$ -mm filters can be placed on the upper surface of the field lens. (If the thickness is more than 2.5 mm, only one filter is allowed.)

## 5

## Troubleshooting

Misuse of this product may adversely affect performance, even if this product is properly functional. If any of the following problems occurs, be sure to check the following table for possible causes before requesting service.

If you detect problems that are not listed below or the problem still persists after measures are taken, unplug the power cord of the device and contact your nearest Nikon representative.

## 5.1 Optical

## Dirt or dust in the viewfield

Problems	Check Items
The specimen is dirty if dirt or dust moves when the specimen is moved on the stage.	<input type="checkbox"/> Clean the specimen.
Objective is dirty if dirt or dust disappears when the objective is switched.	<input type="checkbox"/> Clean the objective.
The eyepiece is dirty if dirt or dust does not move when the stage or objective is moved.	<input type="checkbox"/> Clean the eyepiece.
Aperture diaphragm is closed too far.	<input type="checkbox"/> Open properly. (→ Section 3.5)

## Poor image quality, low resolution, or contrast too low or too high

Problems	Check Items
No coverglass is attached to the slide.	<input type="checkbox"/> Attach a coverglass 0.17 mm thick.
Coverglass is too thick or too thin.	<input type="checkbox"/> Use a coverglass of the specified thickness (0.17 mm).
Slide is upside down.	<input type="checkbox"/> Turn over the slide so that the cover glass faces up.
Lenses or specimen are dirty or dusty.	<input type="checkbox"/> Check dirty or dusty parts, and then clean them as appropriate by reference to the problems described in "Dirt or dust in the viewfield". (→ Section 6.1)
Aperture diaphragm is opened or closed too far.	<input type="checkbox"/> Close or open properly. (→ Section 3.5)
No immersion oil is applied to the tip of the oil-immersion objective.	<input type="checkbox"/> Apply Nikon immersion oil to the objective. (→ Section 3.6)
Nikon immersion oil is not used for oil-immersion observation.	<input type="checkbox"/> Apply Nikon immersion oil to the objective. (→ Section 3.6)
Air bubbles in immersion oil.	<input type="checkbox"/> Remove bubbles. (→ Section 3.6)
Immersion oil adhering to the tip of the dry type objective (especially 40x objective).	<input type="checkbox"/> Clean the lens. (→ Sections 3.6, 6.1)
No coverglass is attached to the slide for observation.	<input type="checkbox"/> Attach a coverglass. 160/-: Can be used for both specimens with and without a coverglass.

**Image too bright**

Problems	Check Items
Setting of the brightness control knob	<input type="checkbox"/> Adjust the brightness by rotating the brightness control knob.

**Insufficient brightness**

Problems	Check Items
Setting of the brightness control knob	<input type="checkbox"/> Adjust the brightness by rotating the brightness control knob.
Aperture diaphragm is closed too far.	<input type="checkbox"/> Open properly. (→ Section 3.5)

**Darkness at the periphery, no viewfield seen, or uneven viewfield brightness**

Problems	Check Items
Revolving nosepiece is not in click-stop position. (the objective is not centered in the optical path).	<input type="checkbox"/> Revolve to click-stop position (swing the objective correctly into the optical path).
Dirt or dust on the lens (condenser, objective, field lens, eyepiece), specimen	<input type="checkbox"/> Clean dirt and dust. (→ Section 6.1)

**Image dark on one side  
(One side of the viewfield (up, down, right, or left) is not focused.)**

Problems	Check Items
Revolving nosepiece is not in click-stop position.	<input type="checkbox"/> Revolve to click-stop position.
Specimen rises from stage surface.	<input type="checkbox"/> Stabilize it using the holder.
The thickness of the sample (section) is not uniform and has a tilt.	<input type="checkbox"/> Prepare a sample (section) of uniform thickness.

**Image shifts during focus  
(Becomes asymmetrically defocused when moving the focal point.)**

Problems	Check Items
Revolving nosepiece is not in click-stop position.	<input type="checkbox"/> Revolve to click-stop position.
Specimen rises from stage surface.	<input type="checkbox"/> Stabilize it using the holder.
The thickness of the sample (section) is not uniform and has a tilt.	<input type="checkbox"/> Prepare a sample (section) of uniform thickness.

## 5.2 Operational

### Image cannot be focused with high magnification objectives.

Problems	Check Items
Slide is upside down.	<input type="checkbox"/> Turn over the slide so that the coverglass faces up.
Coverglass is too thick.	<input type="checkbox"/> Use a coverglass of the specified thickness (0.17 mm).

### High magnification objective contacts slide when changed over from low magnification.

Problems	Check Items
Slide is upside down.	<input type="checkbox"/> Turn over the slide so that the coverglass faces up.
Coverglass is too thick.	<input type="checkbox"/> Use a coverglass of the specified thickness (0.17 mm).
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust. (→ Section 3.4)
Objective is not attached correctly.	<input type="checkbox"/> Screw the objective all the way in. (→ Section 4.4)

### Difference in focal point too large when switching from one objective to another

Problems	Check Items
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust. (→ Section 3.4)
Objective is not attached correctly.	<input type="checkbox"/> Screw the objective all the way in. (→ Section 4.4)

### Image cannot be focused.

Problems	Check Items
The stage vertical movement stopper is lower than the focal position.	<input type="checkbox"/> Turn the stage stopper counterclockwise until it reaches the limit. (→ Section 3.8)

### Binocular images are not integrated.

Problems	Check Items
Interpupillary distance is not adjusted correctly.	<input type="checkbox"/> Adjust. (→ Section 2.6)
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust. (→ Section 3.4)

### Excessive eye fatigue during observation

Problems	Check Items
Interpupillary distance is not adjusted correctly.	<input type="checkbox"/> Adjust. (→ Section 2.6)
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust. (→ Section 3.4)
Inadequate brightness or illumination.	<input type="checkbox"/> Adjust brightness using the brightness control knob.

## 5.3 Electrical

### Illumination does not turn on when the light ON/OFF switch is turned on.

Problems	Check Items
The AC adapter is not connected to the microscope main body.	<input type="checkbox"/> Insert the DC cord of the AC adapter into the DC input connector. (→ Section 2.2)
No electrical power.	<input type="checkbox"/> Insert the power cord connected with the AC adapter into an electrical outlet. (→ Section 2.2)

### Flickering or unstable lamp brightness

Problems	Check Items
The AC adapter or power cord is not correctly connected.	<input type="checkbox"/> Connect correctly. (→ Sections 2.2, 4.6)

# 6

## Care and Maintenance

### 6.1 Cleaning the Lenses

- Dust is best removed using a soft brush or gauze.
- More persistent dirt, such as fingerprints, grease and oil, may be removed with clean lens tissue (or soft cotton, gauze) lightly moistened with absolute alcohol (ethyl alcohol or methyl alcohol).
- To clean immersion oil off the oil-immersion type objective, use lens tissue, soft cotton or gauze lightly moistened with petroleum benzine. If petroleum benzine is not available, use methyl alcohol. In this case, you need to wipe 3 to 4 times because the detergency of the methyl alcohol is somewhat weak.
- Absolute alcohol and petroleum benzine are quite inflammable. Take great care when handling them and when turning the power switch on and off. Be very careful with fire.
- When using absolute alcohol and petroleum benzene, follow the instructions provided by the manufacturer.

### 6.2 Cleaning the Microscope

- We recommend that you use a silicon cloth to clean the microscope.
- For persistent dirt, dampen a piece of gauze with neutral detergent and wipe lightly.
- Do not use organic solvent (alcohol, ether, or thinner). Doing so may result in discoloration of the coating and plastic parts or removal of the printed text.

### 6.3 Disinfecting the Microscope

- We recommend that you use 70% medical alcohol for normal disinfection of the microscope.
- In case of spillage of a sample to the microscope, determine whether the sample is hazardous. If the sample is hazardous, follow the standard procedure of your laboratory.
- Using organic solvent may result in discoloration of the plastic parts.

### 6.4 Storing the Microscope

- When the microscope is not in use, put a vinyl dust cover over the product to protect it from dust, and store it in a dry place where mold is not likely to form.
- Before covering the microscope with the vinyl cover, unplug the power cord.
- We especially recommend that the objectives and eyepieces be kept in a container (such as a desiccator) with desiccant in it.

### 6.5 Periodical Inspections (Charged)

- To maintain the performance of the microscope, periodical inspections and maintenance are recommended.
- For details, contact your nearest Nikon representative.

## 7

## Specifications

## 7.1 Microscopy (Principles)

Use the objective and eyepieces of the microscope to magnify minute cells and tissue optically, manipulate the lever and knobs of the microscope unit to adjust the focus or move the observation point. Then, observe the sample affixed to the slide.

■ **Intended use of the product (for medical care)**

This microscope is intended for use in microscopic experiment and diagnostics of cells and tissues at hospitals or by doctors in private practice in the field of pathology, anatomy, and cytology.

The microscopy with diasopic illumination is used to observe a sample fixed on the slide (cells and tissue) as the specimen.

The product is classified as an in-vitro diagnostic medical device. IVD

This product is not intended for use for measurement.

The scale on the focus knob and stage is an indicator to reproduce the position and does not guarantee the value of the thickness or length of a sample measured using this scale.

■ **Intended user**

It is intended for the medical professional and those who work on experimentations in the field of pathology and cytology.

## 7.2 Performance Properties

<b>Model Name:</b>	Main Body: ECLIPSE Ei R
<b>Optical System:</b>	CFI45 (Infinity-corrected CF optical system) Second objective focal length $f = 200$ mm Built-in diasopic illumination system (Simplified Kohler's illumination system)
<b>Focusing Mechanism:</b>	Fine focus knob graduation: 2 $\mu$ m/graduation Fine focus knob travel: 0.2 mm up or down / revolution Coarse focus knob travel: 37.7 mm up or down / revolution Stage vertical movable range: From the focal plane 2 mm upward and 13 mm downward
<b>Stage:</b>	Stroke: Right and left axis: 76 mm Front and back axis: 30 mm
<b>Revolving Nosepiece:</b>	4-hole fixed type
<b>Illumination Light Source:</b>	LED

## 7.3 Physical Properties

<b>Model Name:</b>	Main Body: ECLIPSE Ei R
<b>External dimensions and Mass (reference)*:</b>	<p>External dimensions: 198 (W) x 377 (H) x 259 (D) mm (Excluding protrusions)</p> <p>Mass (reference): Approximately 5.2 kg</p> <p>* For a set of Ei R and EC-T-TB Binocular eyepiece tube</p> <p>There are some differences in dimensions and mass depending on the details of the set.</p>
<b>Electrical shock protection class</b>	Class III
<b>Environmental conditions</b>	<p>Operation: Temperature: 0 to +40°C Humidity: 60% RH Max. (at +40°C) (no condensation)</p> <p>Storage and Transport: Temperature: -20 to +60°C Humidity: 90% RH Max. (no condensation)</p> <p>Altitude: 2,000m Max.</p> <p>Pollution: Degree 2</p> <p>Overvoltage category: Category II</p> <p>Indoor use only</p>
<b>Safety standards</b>	<ul style="list-style-type: none"> <li>• CE Marking <ul style="list-style-type: none"> <li>• IVD Directive</li> </ul> <p>This equipment complies with the emission and immunity requirements of IEC/EN 61326-2-6.</p> </li> <li>• Low Voltage Directive</li> <li>• EMC Directive</li> <li>• C-UL-US Listed</li> <li>• FCC Part 15 Subpart B Class A</li> </ul> <p>Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.</p> <p>This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.</p> <ul style="list-style-type: none"> <li>• CAN ICES-003(A) / NMB-003(A)</li> <li>• Australian EMI (AS/NZS CISPR11)</li> </ul>


### ■ AC adapter

<b>Name</b>	AC adapter
<b>Manufacturer's name or trade mark, commercial registration number and address</b>	Adapter Technology Co., Ltd. 6F-9 No. 258 Liancheng Rd Zhonghe District New Taipei 235 TAIWAN
<b>Model identifier</b>	ATS018T-P050
<b>Input voltage, input current</b>	Universal 100-240 VAC, single phase, 0.48A Max.
<b>Input AC frequency</b>	50-60 Hz
<b>Output voltage</b>	5.0 VDC
<b>Output current</b>	3.0 A
<b>Output power</b>	15.0 W
<b>Average active efficiency</b>	82.6%
<b>Efficiency at low load (10%)</b>	79.6%
<b>No-load power consumption</b>	>0.04 W
<b>External dimensions</b>	50 (W) x 33 (H) x 100 (D) mm
<b>Mass (reference)</b>	Approx. 170 g (without the power cord)
<b>Electrical shock protection class</b>	Class I
<b>Environmental conditions</b>	<p>Operation:                      Temperature: 0 to +40°C   Humidity: 20 to 80% RH (no condensation)</p> <p>Storage and Transport:      Temperature: -20 to +80°C   Humidity: 10 to 90% RH (no condensation)</p> <p>Altitude:                         2,000m Max.</p> <p>Pollution:                       Degree 2</p> <p>Overvoltage category:        Category II</p> <p>Indoor use only</p>
<b>Safety standards</b>	CE marking, GS mark, C-UL-US Listed, PSE mark


### ■ Power Cord

<b>When used in Japan:</b>	PSE approved detachable power cord set (3 conductor grounding Type VCTF 3x0.75mm <sup>2</sup> , 3m long maximum, rated at 125 VAC minimum)
<b>When used in a 100-120 V region, outside Japan:</b>	UL listed detachable power cord set (3 conductor grounding Type SVT, NO.18 AWG, 3m long maximum, rated at 125 VAC minimum)
<b>When used in a 220-240 V region:</b>	EU/EN standard approved detachable power cord set (3 conductor grounding Type H05VV-F 1 mm <sup>2</sup> , 3m long maximum, rated at 250 VAC minimum)

■ External dimensions

**Binocular eyepiece tube set**

**Trinocular eyepiece tube set**

**En** Anti-fungal measures

- In order to protect the inside of this product from fungi, the objectives, eyepieces and tubes are treated with germicide.
- The active substance used is IPBC.
- Active ingredients are classified as causing skin sensitization.

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**En Symbol for separate collection applicable in European countries**

This symbol indicates that this product is to be collected separately. The following apply only to users in European countries.

- This product is designated for separate collection at an appropriate collection point. Do not dispose of as household waste.
- For more information, contact the retailer or the local authorities in charge of waste management.

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1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A.
NIKON INSTRUMENTS EUROPE B.V.
NIKON INSTRUMENTS(SHANGHAI)CO.,LTD.
NIKON SINGAPORE PTE LTD
NIKON INSTRUMENTS KOREA CO.,LTD.
NIKON FRANCE S.A.S.